beurer
AUTO 400
Upper Arm Blood Pressure Monitor

For Canada only / Solo para Canada /
Pour le Canada seulement:
Manufacturer / Fabricante / Fabricant :
Beurer GmbH, Söflinger Str. 218,
89077 Ulm, Germany

Distributed by / Distribuido por / Distribue par :
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020, USA
www.shop-beurer.com

ENGLISH

1. IMPORTANT SAFETY NOTES

Signs and symbols

Whenever used, the following signs identify safety and property damage measures and designate a level of hazard or seriousness.

READ THIS ENTIRE MANUAL, THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

	This is the safety alert symbol. It alerts you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
	Indicates a hazardous situation which, if not avoided, could result in death or serious injury.
	Indicates a hazardous situation which, if not avoided, may result in minor or moderate injury
	Addresses practices not related to personal injury, such as product and/or property damage.
	Isolation of applied parts Type BF Galvanically isolated applied part (IF stands for "floating"); meets the requirements for leakage currents for type B
	Observe the instructions Read the instructions before starting work and/or operating devices or machines
	Direct current The device is suitable for use with direct current only
	Product information Note on important information
	Temperature limitation
	Humidity limitation
	Serial number
	Protected against solid foreign objects 12.5 mm in diameter and larger, and against vertically falling drops of water

WARNING

- Do not use the blood pressure monitor on newborns or patients with pre-eclampsia. Consult your doctor before using the blood pressure monitor during pregnancy.
- This device may have difficulty determining the proper blood pressure for pregnant women and for users with irregular heartbeat, diabetes, poor circulation of blood, kidney problems, or for users who have suffered from a stroke.
- The pulse display is not suitable for checking the frequency of heart pacemakers or defibrillators.
- This device is not intended for use by people (including children) with restricted physical, sensory or mental skills or a lack of experience and/ or a lack of knowledge, unless they are supervised by a person who has responsibility for their safety or they receive instructions from this person on how to use the device. Supervise children around the device to ensure they do not play with it.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access intravascular therapy, or an arteriovenous (AV) shunt.
- Do not use the device at the same time as other medical electrical devices (ME equipment). This could lead to a malfunction of the device and/or an inaccurate measurement.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions remove the cuff from the arm.

CAUTION

- Place the cuff on your upper arm only. Do not place the cuff on other parts of the body.
- Please note that when inflating, the functions of the limb in question may be impaired.
- Do not perform measurements more frequently than necessary. The resulting restriction of the blood flow may cause injury.

NOTICE:

- Do not use the device outside of the specified storage and operating conditions. This could lead to incorrect measurements.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.
- Using the blood pressure monitor outside your home environment or whilst on the move (e.g. whilst travelling in a car, ambulance or helicopter, or whilst undertaking physical activity such as playing sport) can influence the measurement accuracy and cause incorrect measurements.
- The blood pressure monitor is made from precision and electronic components. The accuracy of the measured values and service life of the device depend on its careful handling:
 - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
 - Do not drop the device.
 - Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
 - Only use the cuff included with the delivery or original replacement cuffs. Otherwise incorrect measured values will be recorded.
- Bring the device to room temperature before taking measurements. If the measuring device has been stored near to the maximum or minimum

storage and transport temperature and is brought into an environment with a temperature of 68 °F (20 °C), it is recommended that you wait 2 hours before using the measuring device.

- Never immerse or spill water or any other liquid onto the monitor or any components, otherwise liquid will enter it and cause damage.
- Never attempt to repair, open and/or disassemble the unit (including arm cuff and optional accessory) or adjust it yourself. This may damage the unit and impair the functions. If you need to have the unit repaired, please contact our customer service. Please see warranty for service contact. Before submitting any complaint, first check the batteries and replace them if necessary.
- Only use the device on people who have the specified upper arm measurement for the device.
- Avoid any mechanical restriction, compression or bending of the cuff line.
- To conserve the batteries, the blood pressure monitor switches off automatically if no buttons are pressed for 1 minute.
- Do not expose the device to static electricity. Always ensure that no static electricity is coming off of you before you operate the device.

Battery Handling Safety Precautions

- Use only the size and type of batteries specified.
- Be sure to follow the correct polarity when installing the batteries. Reversed batteries may cause damage to the device.
- Do not mix different types of batteries together (e.g. Alkaline and Carbon-zinc or rechargeable batteries) or old batteries with fresh ones. Always replace batteries as a simultaneous set.
- If the batteries in the device are depleted or the device will not be used for a long period of time, remove the batteries to prevent damage or injury from possible battery leakage.
- Do not try to recharge batteries not intended to be recharged; they can overheat and rupture (follow battery manufacturer's directions.)
- Do not dispose of batteries in fire, batteries may explode or leak.
- Protect batteries from excessive heat.
- Clean the battery contacts and also those of the device prior to battery installation.
- Remove discharged batteries from the product and dispose/recycle in compliance with all applicable laws.
- Keep batteries away from children and pets. Batteries may be harmful if swallowed. Should a child or pet swallow a battery, seek medical assistance immediately.
- If your skin or eyes come into contact with battery fluid, rinse the affected areas with water and seek medical assistance.
- Do not disassemble, open or crush the batteries.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can impair the performance of the device.


2. INTENDED USE

This blood pressure monitor performs non-invasive measurement and monitoring of arterial blood pressure and heart rate on adult patients.

3. PACKAGE CONTENTS


- 1 x Upper Arm Blood Pressure Monitor
- 1 x Universal cuff (8.7 in - 16.5 in / 22 cm - 42 cm)
- 4 x 1.5 V AAA batteries
- 1 x Quick Start Guide
- 1 x Instruction manual

4. DEVICE DESCRIPTION

- Universal cuff
- Air Tube
- Air Tube Connector
- Connection for Air Tube connector
- Risk Indicator
- LCD Display
- MEMORY Button
- START/STOP Button
- > Button
- Battery Compartment Cover

LCD Display Information

- Memory space number
- Time and Date
- Systolic Pressure
- Diastolic Pressure
- Pulse value determined
- Irregular Heartbeat Indicator
- Low Battery Icon
- User memory
- Risk indicator
- Cuff position indicator
- Memory display: average value

5. PREPARING FOR MEASUREMENT

Inserting/changing batteries:

- Remove the Battery Compartment Cover from the rear of the device.
- Insert four "AAA" alkaline batteries (included) according to the polarity markings inside the compartment. **NOTE:** Do not use rechargeable batteries.
- Reattach the Battery Compartment Cover.
- All LCD Display elements will briefly appear and "12h" will flash. Set the date and time as described below. Any saved measurements will be retained.

If the battery replacement symbol is permanently displayed, you can no longer perform any measurements and must replace all batteries. Once the batteries have been removed from the device, the date and time must be set again. Any saved measurements are retained.

Setting the hour format, date, and time:

If date and time are not set, stored measurements will not have a valid timestamp associated with them.

The settings menu can be accessed in two different ways:

Before first use and after every battery change:

When you insert batteries into the device, you will automatically arrive at the relevant menu.

If batteries have already been inserted:

With the device switched-on press and hold the **MEMORY** Button for 3 seconds.

- The time format will flash. Select an hour format (12- or 24-hour display) with the > button and press the **MEMORY** Button to confirm.
- The year flashes on the LCD Display. Set the year with the > button and press the **MEMORY** Button to confirm.
- The month flashes on the LCD Display. Set the month with the > button and press the **MEMORY** Button to confirm.
- The day flashes on the LCD Display. Set the day with the > button and press the **MEMORY** Button to confirm.
- NOTE:** If the hour format is set to 24h, the day/month display sequence is reversed.
- The hour flashes on the LCD Display. Set the hour with the > button and press the **MEMORY** Button to confirm.
- The minutes flash on the LCD Display. Set the minutes with the > button and press the **MEMORY** Button to confirm.

6. MEASURING BLOOD PRESSURE

The blood pressure monitor is suitable for home use and needs to be at room temperature when measuring. Measurements can be performed on either arm.

Positioning the Cuff:

- Place the cuff on a bare upper arm. Ensure blood circulation is not restricted by tight clothing or other objects. The cuff must be placed on the upper arm so that the bottom edge is positioned 0.8 in - 1.2 in above the elbow and over the artery, with the Air Tube in line with the center of the palm.
- Fasten the cuff, making sure it is not too tight around your arm. You should be able to fit two fingers under the cuff after it is in place. The cuff is suitable for you if the index mark is within the OK range on the cuff when it is in place.

- Insert the Air Tube Connector into the Air Tube Port.


Blood pressure can vary between the right and left arm, which may mean that the measured blood pressure values are different. Always perform the measurement on the same arm. If the values between the two arms are significantly different, please consult your doctor to determine which arm should be used for the measurement.

Cuff Position Indicator

This device features a cuff position indicator to check the correct positioning of the cuff on your upper arm during measurement.

When the cuff is correctly positioned, the symbol will appear in the LCD Display during measurement and will also display along with measurement results. If the cuff is positioned incorrectly, the measurement is cancelled and an error message appears in the LCD Display. If this happens, refer to the "Troubleshooting Guide" section for solutions.

Correct Posture:

- Rest for at least five minutes before each measurement.
- Perform the measurement either sitting or lying down. Always make sure that the cuff is level with your heart.
- To perform a measurement, sit comfortably with your arms and back supported. Do not cross your legs and place your feet flat on the floor.
- Keep still during the measurement and do not talk.
- Relax your arm and the palm of your hand.

Performing a Measurement

- With the cuff attached properly and in a correct posture, start the device using the **START/STOP** Button . Press the **START/STOP** Button again; the device will then automatically begin measuring by inflating the cuff.
- The measurement is taken while the cuff is being inflated. As soon as a pulse has been detected, the pulse symbol will flash in the rhythm of your pulse. After no pulse is detected anymore, the cuff air pressure will be released.
- The measurements for systolic pressure, diastolic pressure and pulse will be displayed. You can stop the measurement at any time by pressing the **START/STOP** button .
- The measurement will be saved automatically.

- "Er..." will appear if the measurement was not performed properly. Refer to the "Troubleshooting Guide" section for details about specific error message numbers and then repeat the measurement.

Note: The device switches off automatically after 60 seconds if no button has been pressed.

Wait for at least five minutes before taking another measurement!

7. EVALUATING RESULTS

Irregular Heartbeat:

This device can detect irregular heartbeat patterns during measurement and will indicate them by displaying the Irregular Heartbeat Icon . This could be an indicator for arrhythmia, a condition where the heart rhythm is abnormal. These symptoms (omitted or premature heartbeats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress, or lack of sleep. This function is not designed for diagnosing or treating an arrhythmic disorder. Only your doctor can determine if you do have arrhythmia.

Repeat the measurement if the Irregular Heartbeat Icon displays after a measurement. Rest for at least five minutes between measurements and do not talk or move during measurement. If the icon appears often, contact your doctor, since self-diagnosis and treatment based solely on the test results may be dangerous. It is vital to follow your doctor's instructions.

Risk Classification

Measurements can be classified and assessed by the table below. However, these standard values are only a general guideline, since individual blood pressure varies in different people and different age groups, etc.

Hypertension Category	Systolic (in mmHg)	Diastolic (in mmHg)	LED Risk Indicator
Grade 3: Severe hypertension	≥ 180	≥ 110	Red
Grade 2: Moderate hypertension	160–179	100–109	Orange
Grade 1: Mild hypertension	140–159	90–99	Yellow
High Normal	130–139	85–89	Green
Normal	120–129	80–84	Green
Optimal	< 120	< 80	Green

Source: WHO, 1999

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

The classification on the display and the scale on the unit show which category the recorded blood pressure values fall into. If the systolic and diastolic values fall into two different categories, the Risk Classification Indicator on the device always shows the higher category.

8. DISPLAYING AND DELETING MEASUREMENTS

User Memory:

Every successful measurement is stored along with its date and time. If there are more than 60 measurements, the oldest measurement is deleted.

Average Values:

- With the device switched on, press the **MEMORY** button.
- The average value of all measurements will appear on the display.
- Press the > button to display the average value of the last 7 days for the morning measurement. (Morning: 5:00 am–9:00 am, display). Press the > button again to display the average value of the last 7 days for the evening measurement. (Evening: 6 pm–8 pm, display). By pressing the > button again, the last individual measurement will be displayed with the date and time. Switch through all the saved measurements by pressing the > button again.
- Press the **START/STOP** button for 2 seconds to switch off the device.

Clearing Memory:

- To delete the memory, press the **MEMORY** button when the device is switched on.
- The average value of all measurements will appear on the display.
- Keep the > button and **MEMORY** button pressed simultaneously for 5 seconds.
- After successful deletion, "CL 01" will appear in the display. The data has been deleted.

9. CARE AND MAINTENANCE

- Clean the device and cuff with a slightly moistened cloth.
- Do not use detergents or solvents.
- Never immerse the device in water or other liquids.
- When storing the device, make sure no heavy objects are placed on top of it. Remove the batteries. Ensure the Cuff line does not have any sharp kinks.

Disposal

Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

10. TROUBLESHOOTING GUIDE

In case of errors, "Er..." will appear in the LCD Display, where _ is a number 1-6. Error messages may appear if:

Trouble shooting	Problem	Solution
"Er1"	The pulse cannot be measured correctly.	Repeat the measurement process. Be sure not to move or talk.
"Er2"	You have been talking or moving during the measurement process.	
"Er3"	The cuff was attached too tightly or too loosely.	
"Er4"	An error occurred during the measurement process.	
"Er5"	The pump-up pressure is higher than 300mmHg or the blood pressure measurement falls outside of the measuring range.	Repeat the measurement process. Be sure not to move or talk.
"Er6"	A device error has occurred.	Repeat the measurement process. Be sure not to move or talk.
L0	The batteries are almost dead.	Reinsert the batteries or replace them with new batteries.

11. TECHNICAL SPECIFICATIONS

Type	AUTO 400
Model	BM 72
Measurement Method	Oscillometric, non-invasive blood pressure measurement on the upper arm
Measurement Range	Measuring range 30-280 mmHg Systolic 50 – 280 mmHg Diastolic 30 – 200 mmHg Cuff pressure 0 – 300 mmHg Pulse 40 – 200 beats per minute
Display Accuracy	Systolic ± 3 mmHg, Diastolic ± 3 mmHg, Pulse ± 5% of the displayed value
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: systolic 8 mmHg/diastolic 8 mmHg
Memory	60 memory spaces
Dimensions	L 4.92 in x W 3.74 in x H 1.65 in (L 125 mm x W 95 mm x H 42 mm)
Weight	0.73 lb (335 g) with cuff, without batteries
Cuff Size	8.7 in - 16.5 in (22 - 42 cm)
Operating Range	50°F - 104°F (10 °C - 40 °C), <90% relative humidity (non-condensing), 800-1050 hPa ambient pressure
Storage Range	-4 °F - 131 °F (-20 °C - 55 °C), ≤ 90% relative humidity (non-condensing), 800-1050 hPa ambient pressure
Power Supply	4 x 1.5V === AAA batteries
Battery Life	Approx. 350 measurements

The serial number is located on the device or in the battery compartment. Technical information is subject to change without notification to allow for updates.

- This unit is in line with IEC 80601-2-30 Safety and Performance, IEC 60601-1 Electrical Safety and IEC 60601-1-2 Electromagnetic Compatibility.
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life.

12. FCC COMPLIANCE INFORMATION

Upper Arm Blood Pressure Monitor BM 72

Responsible Party – U.S. Contact Information

Beurer North America LP

1 Oakwood Boulevard, Suite 255

Hollywood, FL 33020

United States

1-800-536-0366

info@beurer.com

FCC Compliance Statement

This device complies with Part 15 of the FCC Rules.Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

13. WARRANTY

Limited Lifetime Warranty For Original Purchaser

Your Beurer AUTO 400 Upper Arm Blood Pressure Monitor, model BM 72, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the Beurer AUTO 400 Upper Arm Blood Pressure Monitor, model BM 72, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.

For warranty service contact our customer service department at 1-800-536-0366 or at info@beurer.com to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessory; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized repairs or modifications; improper use of electrical/power supply; old worn batteries; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Beurer. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of Beurer under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY.

Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL BEURER BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you. Beurer does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate Beurer in any way beyond the terms set forth herein. This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resaled products, including but not limited to sale of such products on Internet auction sites and/or products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of Beurer.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state. For more information regarding our product line in the USA, please visit: www.shop-beurer.com

For Canada only:

Manufacturer:

Beurer GmbH

Söflinger Str. 218

89077 Ulm Germany

Distributed by:

Beurer North America LP

1 Oakwood Boulevard, Suite 255

Hollywood, FL 33020

USA

www.shop-beurer.com

Questions or comments? Call our US-based customer service toll free at 1-800-536-0366.

Assembled in China

EMC Guidance

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments

Warning: Don't use near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

- All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.
- Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions

Emissions test	Compliance
RF emissions, CISPR 11	Group 1,
RF emissions, CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions, IEC 61000-3-3	Not applicable

	Test fre- quency (MHz)	Band (MHz)	Service	Modulation	Maximum Power (W)	Distance (m)	Test Level (V/m)	Compliance level (V/m)
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equip- ment)	385	380-390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27	27
	450	430-470	GMRS 460, FRS 460	FM ± 5 kHz de- viation 1 kHz sine	2	0,3	28	28
	710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9	9
	745							
	780							
	810							
	870							
	930	900-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28	28
	1720							
	1845							
	1970							
	2450	2400- 2570	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28	28
	5240 5500 5785	5100- 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9	9

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