beurer

FC 76

ENGLISH **Microdermabrasion**
Instruction for Use 3

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Questions or Comments? Call our US-based customer service toll free at 1-800-536-0366.

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READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT

Keep this manual in a safe location for future reference

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ENGLISH

Contents

1. IMPORTANT SAFETY NOTES4
2. Getting to Know your Instrument6
3. Package Contents6
4. Parts and Controls.....7
5. Before Use.....7
6. To Use7
7. Care, Maintenance and Disposal ..10
8. Replacement Parts.....10
9. Troubleshooting Guide10
10. Technical Specifications.....11
11. FCC Compliance Information.....11
12. Warranty12

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE
If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact Beurer customer service. See warranty for service contact.

Dear Valued Customer,
Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for application in the areas of heating, weight management, blood pressure, body temperature, pulse, gentle therapy, massage, and beauty.

Register your product today to receive further benefits:
- added protection for your purchase
- health and wellness tips
- new product information
- continual outstanding customer service and support.

Visit www.registerbeurer.com or contact customer service at 1-800-536-0366.
Please read these instructions for use carefully and keep them for later use.
Make them accessible to other users and observe the information they contain.

With kind regards,
Your Beurer team

1. IMPORTANT SAFETY NOTES

Signs and symbols

Whenever used, the following signs identify safety and property damage messages and designate a level of hazard or seriousness.

READ THIS ENTIRE MANUAL, THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

	This is the safety alert symbol. It alerts you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.		Polarity
WARNING	Indicates a potentially hazardous situation, which, if not avoided, could result in death or serious injury.		Only use indoors
NOTICE	Addresses practices not related to personal injury, such as product and/or property damage.		**DANGER** This device must not be used near to water or in water (e.g. in a wash basin, or in the shower or bath) – risk of electric shock!
CAUTION	Indicates a potentially hazardous situation, which, if not avoided, may result in minor/moderate injury.	**IP20**	Protected against solid foreign objects 12.5 mm in diameter and larger
	**Product information** Note on important information		

⚠ WARNING – To reduce the risk of fire, electric shock, or serious personal injury:

- Danger of suffocation by packaging material. Keep packaging material away from children.
- Before use, inspect the device and accessories for visible damage. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.
- The device is intended solely for indoor, non-commercial, non-industrial household use on facial skin (apart from the eyelids and lips) of human adults. Do not use this device for any other purpose.
- This device is not suitable for persons under the age of 18.
- Ensure children do not play with the device.
- Cleaning and user maintenance must not be performed by children.
- If the power cord of this device is damaged, discard it immediately and replace it with a new cable.
- Do not open or repair the device yourself; there are no user serviceable parts inside. Opening the device will void the warranty.
- Do NOT use this device:
 - on cracked skin
 - on open wounds
 - in the eye area or on the mouth
 - on any part of the body other than the face (apart from eye area and mouth)
 - during pregnancy
 - if you have had facial cosmetic surgery
 - if you suffer from any type of heart disease
 - if you have a skin disease (e.g. herpes, vascular diseases or skin growths)
 - if you have skin cancer
 - if you suffer from skin irritation
 - immediately after sunbathing or if you have sunburn
 - if you are taking medication that contains steroids
 - if you are taking blood thinners
 - if you have multiple sclerosis, diabetes, or an autoimmune disorder
- Keep the power cord and power adapter away from water, steam, and other liquids.
- Only use the device in dry rooms (e.g. never in the bath, sauna).
- Do not drop the device or subject it to any impacts.
- Only use the device with dry hands.
- Do not clamp or pinch the power cord.
- Never reach for a device that has been submerged in water. Disconnect the power adapter immediately.
- Only use this device once per week.

- Immediately after use, the skin in the treated area may be red and tender; this is normal. Reddening should subside after about 30 minutes.

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can impair the performance of the device.


2. Getting to Know your Instrument

This device enables professional-quality facial peels with microdermabrasion technology. It uses a high-quality sapphire coating to gently remove dead skin cells as well as state-of-the-art vacuum technology that stimulates skin surface renewal and stimulates skin surface renewal for tighter and cleaner looking skin.
The device is only intended for private use, and not for medical or commercial purposes. Use the device for the purpose it was developed for only and in the manner specified in these instructions for use. Any form of improper use can be dangerous. The manufacturer is not liable for damage resulting from improper or incorrect use. Do not use the device if one or several of the following warnings and safety notes applies to you. If you are unsure whether the device is suitable for you, consult your doctor.

3. Package Contents

1 x Microdermabrasion FC 76
1 x Precision attachment
1 x Coarse attachment
1 x Fine attachment
20 x Single-use filters
1 x AC power adapter
1 x Instruction manual

4. Parts and Controls


1. Precision attachment with fine sapphire coating - For targeted treatment of small areas around the face
2. Fine attachment with fine sapphire coating - For sensitive skin
3. Coarse attachment with coarse sapphire coating - For a more in-depth treatment
4. Hand held device
5. Power Switch

5. Before Use

About 24 hours before using the device for the first time, perform a test with the device to ensure your skin is not too sensitive for treatment. Test the device at the lowest intensity level on a small, inconspicuous area of your neck. If, after 24 hours, the tested area is not red or painful, you may continue using the device. Otherwise, discontinue use immediately.


6. To Use

1. Connect the power cord to the power jack on the bottom of the Base and to a standard electrical outlet.


2. Make sure your face is clean, dry, and free of makeup, creams, and powders.


3. Place a filter on the Base and replace the Filter Cover.


4. Select an Attachment and attach it to the Base.

5. Slide the Power Switch to the desired intensity level (Low or High).

6. Move the device around your face, making sure to treat each area only once. For details on treatment techniques, see the "Application Areas" subsection below. Also make sure not to apply excess pressure during treatment, as using this device should be comfortable and not painful. Even so, some mild reddening of the treatment area is normal.
7. When treatment is finished, slide the Power Switch to the Off position and unplug the power cord.
8. Change the Filter after each use by removing the Nozzle and Filter Cover and using tweezers to remove the filter. Replace with a new Filter and reattach the Filter Cover and Nozzle.


9. For one day before and after each use, avoid activities that cause significant perspiration, such as saunas, tanning beds, and sun exposure. Gently wash and dry your face after each treatment session. Most people should use this device only once per week, but it can be used up to twice a week on non-sensitive problem skin if allowing at least a 48-hour break in between treatment sessions.

Application Areas

Ensure the treatment area skin is taut during treatment. To accomplish this, simply pull the skin gently to one side or downward with your free hand. Do not use the device on one area for a long period.

Forehead:

Move the device slowly and evenly out from the center of your forehead.

Temples:

Move the device outward toward your hairline. Repeat on other side. Keep device at least ½ inch (about 1 cm) away from your eyes.

Nose:

Move the device down the center of your nose starting at the top. Repeat for the left and right sides.

Cheeks:

Move the device outward from the nose across your cheeks.

Mouth Area:

Move the device along your upper lip from the middle to the outside. Then move the device up and outwards from your bottom lip along the jawline. Make sure that you do not touch your lips.

7. Care, Maintenance and Disposal

- Never submerge the device in water or other liquids.
- Disinfect the attachments with a standard disinfectant after each use.
- Change the filter after each use.
- Use a dry cloth to clean the device of any dirt. Do not use any abrasive cleaning products.
- Protect the device from impacts, moisture, dust, chemicals, large temperature fluctuations and nearby sources of heat (ovens, heaters).

Disposal

Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

8. Replacement Parts

To order replacement parts, contact Customer Service and reference the following part numbers.

Part Description	Part Number
10 Replacement Filters	163.660
Precision Attachment	163.767
Coarse Attachment	163.766
Large Attachment	163.765

9. Troubleshooting Guide

Problem	Cause	Solution
Device will not switch on	On/Off switch faulty	Contact Customer Services or your retailer
No power	The AC power adapter is not connected correctly.	Connect the AC power adapter correctly.
No power	The AC power adapter is faulty.	Contact Customer Services or your retailer
The device's suction is incorrect	Filter blocked/pump faulty	Replace the old filter with a new filter
Exfoliation effect from attachments is no longer ensured	Sapphire coating has worn down	Contact Customer Services or your retailer to request new attachments.

10. Technical Specifications

Device:	
Type	FC 76
Power supply	DC 5.0 V, 1.0 A
Dimensions	7.5 x 1.6 x 1.5 inches (190 x 40 x 39 mm)
Weight	4.5 oz (128g)

AC power adapter:	
Modelnumber	LGSPSD050100UL
Input	AC 100-240 V, 50/60 Hz, 0.5 A
Output	DC 5.0 V, 1.0 A, 5.0 W
	Polarity of the DC voltage connection
Supplier	Lygo (Dongguan) Electronic Product Ltd
Safety Standard	UL 1310 Class 2 power unit and CAN/CSA C22.2 No. 223 - Power Supplies with Extra low Voltage class 2 Output

11. FCC Compliance Information

Microdermabrasion FC 76

Responsible Party – U.S. Contact Information
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020, United States
1-800-536-0366, info@beurer.com

FCC Compliance Statement
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
Changes or modifications to the product not expressly approved by the manufacturer could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to

provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency
energy and, if not installed and used in accordance with the instructions, may cause harmful interference to
radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one
or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

12. Warranty

Limited Lifetime Warranty For Original Purchaser

Your Beurer Microdermabrasion, model FC 76, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the Beurer Microdermabrasion, model FC 76, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.

For warranty service contact our customer service department at 1-800-536-0366 or at info@beurer.com to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessory; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized

repairs or modifications; improper use of electrical/power supply; old worn batteries; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Beurer. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of Beurer under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL BEURER BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you. Beurer does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate Beurer in any way beyond the terms set forth herein.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of Beurer.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state. For more information regarding our product line in the USA, please visit: www.beurer.com

Distributed by:
Beurer North America LP
1 Oakwood Boulevard, Suite 255, Hollywood, FL 33020, USA
www.beurer.com

Questions or comments? Call our US-based customer service toll free at 1-800-536-0366.

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