beurer

EM 46

ENGLISH **Digital TENS/EMS Device**
Instructions for use 3

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READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT

Keep this manual in a safe location for future reference

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ENGLISH

CONTENTS

1. IMPORTANT SAFETY INSTRUCTIONS....3
2. Indications for Use....4
3. General warnings....6
4. General Information....9
5. Package Contents....9
6. Parts and Controls....10
7. Initial Setup....10
8. To Use....11
9. Program Overview....12
10. Favorite Program....13
11. Therapy Memory....13
12. Electrode Positioning....14
13. Electric Current Parameters....16
14. Care, Maintenance and Disposal....17
15. Troubleshooting Guide....17
16. Technical Specifications....18
17. Replacement and Consumable Part....18
18. FCC Compliance Information....19
19. Warranty....19

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact Beurer customer service. Please see warranty for service contact.

1. IMPORTANT SAFETY INSTRUCTIONS

Signs and symbols

Whenever used, the following signs identify safety and property damage messages and designate a level of hazard or seriousness.

READ AND UNDERSTAND THIS ENTIRE MANUAL INCLUDING THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

⚠	This is the safety alert symbol. It alerts you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
⚠ **WARNING**	Indicates a hazardous situation which, if not avoided, could result in death or serious injury.

⚠ **CAUTION**	Indicates a hazardous situation which, if not avoided, may result in minor or moderate injury.
NOTICE	Addresses practices not related to personal injury, such as product and/or property damage

	**Product information** Note on important information
	**Observe the instructions** Read the instructions before starting work and/or operating devices or machines
IP22	Device protected against foreign objects ≥ 12.5 mm and against water dripping at an angle
SN	Serial number
	Type BF applied part
	Temperature limit
	Humidity limit
	The device may not be used by people with medical implants (e.g. pacemakers), as this may affect their functionality.

	Magnetic Resonance (MR) unsafe. Keep away from magnetic resonance imaging (MRI) equipment. The device presents a projectile hazard.
	This device can emit effective output values above 10 mA, averaged over every 5-second interval

2. INDICATIONS FOR USE

TENS: the device is designed to be used for temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household and work activities.
EMS: The device is designed to be used for stimulating healthy muscles in order to improve and facilitate muscle performance.

⚠ CONTRAINDICATIONS

Do not use the device near the heart. Never place stimulation electrodes on the front ribcage (where the ribs and breastbone are located), and especially not on the two large pectoral muscles. This can increase the risk of ventricular fibrillation and induce cardiac arrest (see “General warnings” section).

- Do NOT use this device
 - On the skeletal skull structure, or around the mouth, throat or larynx.
 - In the neck/carotid artery area.
 - In the genital area.
 - If you have implanted electrical devices (e.g. a pacemaker).
 - If you use a cardiac demand pacemaker.
 - If you have metal or electrical implants.
 - If you use an insulin pump.
 - If you have a high temperature (e.g. > 102.2°F / 39°C).
 - If you have a known or acute cardiac arrhythmia or disorders of the heart’s conduction system.

- If you suffer from a seizure disorder (e.g. epilepsy).
- If you are pregnant.
- For relieving pain associated with abdominal cramping or dysmenorrhea, or over a menstruating or pregnant uterus.
- If you have cancer.
- After an operation, if strong muscle contractions could affect the healing process.
- If you are connected to a high-frequency surgical device.
- If you have an acute or chronic disease of the gastrointestinal tract.
- In the case of a known allergy to the electrode material.
- On acutely or chronically affected (by injury or irritation) skin (e.g. inflamed skin – whether painful or not, reddened skin, rashes, e.g. due to allergies, burns, bruising, swelling, open and healing wounds and post-operative scars where the healing process could be impaired).

• The TENS device is MR Unsafe. Keep away from magnetic resonance imaging (MRI) equipment. The device presents a projectile hazard

⚠ WARNING

Undesirable side effects and adverse reactions

- Skin irritation and burns beneath the electrodes have been reported with the use of TENS.
- Feeling of pressure at the electrode location
- Slight redness, burning and pain of the skin after treatment
- Paraesthesia
- Discomfort
- Sleepiness
- Vibration of muscles
- Tension
- Headaches
- Increased menstrual bleeding
- Allergic reactions to components

MedWatch is the Food and Drug Administration‘s (FDA) program for reporting serious reactions, product quality problems, therapeutic inequivalence/failure, and product use errors with human medical products, including drugs, biologic products, medical devices, dietary supplements, infant formula, and cosmetics.

If you think you or someone in your family has experienced a serious reaction to a medical product, you are encouraged to take the reporting form to your doctor. Your health care provider can provide clinical information based on your medical record that can help FDA evaluate your report.

However, we understand that for a variety of reasons, you may not wish to have the form filled out by your health care provider, or your health care provider may choose not to complete the form. Your health care provider is not required to report to the FDA. In these situations, you may complete the Online Reporting Form yourself.

You will receive an acknowledgement from FDA when your report is received. Reports are reviewed by FDA staff. You will be personally contacted only if we need additional information.

Submitting Adverse Event Reports to FDA

Use one of the methods below to submit voluntary adverse event reports to the FDA:

Report Online at www.accessdata.fda.gov/scripts/medwatch/index.cfm?action=reporting.home

Consumer Reporting Form FDA 3500B. Follow the instructions on the form to either fax or mail it in for submission. For help filling out the form, see MedWatchLearn. The form is available at https://www.fda.gov/media/85598/download

Call FDA at 1-800-FDA-1088 to report by telephone.

Reporting Form FDA 3500 commonly used by health professionals. The form is available at https://www.fda.gov/media/76299/download

3. GENERAL WARNINGS

- The long-term effects of chronic electrical stimulation are unknown.
- Stimulation should not be carried out:
 - Over the carotid sinus nerves, particularly in patients with a known sensitivity to the carotid sinus reflex.
 - Over the neck or mouth. Severe spasm of the laryngeal and pharyngeal muscles may occur, and the contractions may be strong enough to close the airway or cause difficulty in breathing.
 - Transcerebrally about the brain
 - Over swollen, infected, or inflamed areas or skin eruptions (e.g. phlebitis, thrombophlebitis, varicose veins, ect.
 - Applied over, or in proximity to, cancerous lesions
 - With the electrode surfaces placed on the chest, upper back or across the heart.

⚠ WARNING

The device is not a substitute for medical consultation and treatment. Always consult your doctor first in the event of pain or illness. Before using the device, consult your doctor if any of the following apply to you:

- If you suffer from an acute medical condition, in particular if you suspect or have been diagnosed with high blood pressure, a blood coagulation disorder, propensity to thrombo-embolic conditions or recurrent malignant growths.
- If you have any skin conditions.
- If you have unexplained chronic pain in any part of the body.
- If you suffer from diabetes.
- If you have any sensory impairment that reduces the feeling of pain (e.g. metabolic disorders).
- If you are receiving medical treatment.
- In the event of complaints linked to the stimulation treatment.
- If you suffer from persistently irritated skin due to long-term stimulation at the same electrode site.

ONLY USE THIS DIGITAL TENS/EMS DEVICE:

- On people.
- For the intended purpose and as specified in these instructions for use. Any form of improper use can be dangerous.
- For external use only.
- With the supplied and reorderable genuine replacement parts to ensure safe use. Failure to do so voids the warranty and may result in risks for the user. Only use the replacement parts specified in the accompanying documents.
- This device is NOT intended for use by people with restricted physical, sensory or mental skills or a lack of experience and/or a lack of knowledge, unless they are supervised by a person who has responsibility for their safety or they receive instructions from this person on how to use the device.
- Children must not play with the device.
- The device is only intended for domestic/private use, not for commercial use.

⚠ GENERAL PRECAUTIONS

- Caution should be used for patients with suspected or diagnosed heart problems.

- Caution should be used for patients with suspected or diagnosed epilepsy.
- Caution should be used in the presence of the following:
 - When there is a tendency to hemorrhage following acute trauma or fracture;
 - Following recent surgical procedures when muscle contraction may disrupt the healing process;
 - Over areas of the skin which lack normal sensation.
- Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. The irritation can usually be reduced by using an alternate conductive medium, or alternate electrode placement.
- TENS should be used only with the leads and electrodes recommended for use by the manufacturer.
- Do not use in humid environments (e.g. in the bathroom) or when bathing or showering.
- Do not use after consuming alcohol.
- Always pull gently on electrodes to remove them from the skin to prevent injuries in the unusual case of highly sensitive skin.
- Hold the device away from sources of heat and do not use it in close proximity (approx. 39.37 in / 1 m) to shortwave or microwave devices (e.g. mobile phones), as doing so can result in unpleasant current peaks.
- Do not expose the device to direct sunlight or high temperatures.
- Protect the device from impact, dust, dirt and moisture.
- Never immerse the device in water or other liquids.
- The device is suitable for self-treatment.
- For hygiene reasons, the electrodes may only be used on one person.
- If the device does not work properly, or if you feel unwell or experience pain, stop using it immediately.
- Switch off the device or the respective channel first before removing or moving the electrodes to prevent unintentional stimulation.
- Do not modify electrodes (e.g. by cutting them). Doing so increases the current density, which is potentially hazardous.
- Make sure that the electrodes are in full contact with the skin.
- Using worn-out electrodes can cause skin irritation to occur, because an even distribution of current over the entire surface is no longer guaranteed. For this reason, the electrodes should be replaced if any discolouration occurs.
- Do not use whilst asleep, driving a vehicle or operating machinery.
- Do not use whilst undertaking any activity where an unexpected reaction (e.g. strong muscle contractions even at low intensity) could be dangerous.
- Ensure that no metallic objects (e.g. belt buckles or necklaces) come into contact with the electrodes during stimulation. If you are wearing jewellery or have piercings in the area to be treated (e.g. a navel piercing), these must be removed before using the device; failure to do this could result in spot burns.
- Keep the device away from children to prevent potential hazards.
- Make sure not to confuse the electrode cables and their contacts with your headphones or other devices, and do not connect the electrodes to other devices.
- Do not use this device whilst using other devices that transmit electrical pulses to your body.
- Do not use in the vicinity of highly flammable substances, gases or explosives.
- During the initial few minutes, use the device while sitting or lying down to minimize the risk of accidental injuries as a consequence of isolated cases of vagal responses (feeling of faintness). If you feel faint, switch off the device immediately,

lie down and support the legs in an elevated position (approx. 5–10 minutes).

- Treatment of the skin with moisturising lotions or ointments beforehand is not recommended as this considerably increases electrode wear and may cause uncomfortable current peaks.
- Keep packaging material away from children (risk of suffocation!).
- Store the device in a dry place (use indoors only). To avoid the risk of fire and/or electric shock, protect the device from high humidity and water.
- If the adhesive capability of the adhesive electrodes decreases, please replace them immediately. Do not use the device until you have fitted the new adhesive electrodes. Otherwise the unequal adhesion of the adhesive electrodes may lead to skin injuries. Replace the electrodes with new ones at the latest after having used them 20 times.

Damage

- If the device is damaged, do not use it and contact your retailer or the specified Customer Services address.
- To ensure that the device functions effectively, do not drop it or dismantle it.
- Check the device for signs of wear and tear or damage. If there are such signs of wear and tear or damage or if the device was used improperly, it must be returned to the manufacturer or retailer before further use.
- Switch the device off immediately if it is faulty or not working properly.
- Never attempt to open and/or repair the device yourself. Repairs may only be carried out by Customer Services or authorised retailers. Failure to comply will invalidate the warranty.
- The manufacturer is not liable for damage resulting from improper or incorrect use.

Notes on handling batteries

⚠ WARNING

- **Risk of explosion! Risk of fire!** Failure to comply with the following points can result in personal injury or cause overheating, leakage, venting, breakage, explosion, or fire on the battery.
- This device contains non-rechargeable batteries which must not be charged.
- Do not throw batteries into a fire.
- Never charge, forcibly discharge, heat, disassemble, open, crush, deform, encapsulate, or modify batteries.
- Never short-circuit batteries or battery compartment contacts.
- Protect the batteries from direct sunlight, rain, heat, and water.
- Exposure of batteries to an environment with extremely high temperatures or an extremely low air pressure may result in explosion or leakage of flammable liquids and gases.
- Dispose of defective and discharged batteries immediately and properly (see chapter on disposal).
- Do not use modified or damaged batteries.
- Always select the correct battery type.
- Always insert the batteries correctly, taking into account the polarity (+ / -).
- Never mix batteries of different manufacturers, capacities (new and used), size, or type within a device.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- If fluid from a battery comes into contact with your skin or eyes, wash the affected areas with water and seek medical assistance.

- **Choking hazard!** Keep batteries out of the reach of children. Seek medical attention immediately if swallowed. Swallowing them may cause burns, severe internal injuries, and death.
- Never allow children to replace batteries without adult supervision.

⚠ CAUTION

- Store batteries in a well-ventilated, dry, and cool place in a non-conductive container in which the batteries cannot be short-circuited to each other or by other metal objects.
- Keep batteries clean and dry.
- Keep batteries away from water.
- If the device is not going to be used for a long period of time, remove the batteries from the battery compartment.

Notice

- Do not use rechargeable batteries.

4. GENERAL INFORMATION

TENS (transcutaneous electrical nerve stimulation) is the electrical stimulation of nerves through the skin. TENS is a safe and effective, medication-free method of treating certain types of pain. It is also cleared as a means of simple self-treatment.


In sports and fitness, electrical muscle stimulation (EMS) is used to complement conventional muscle training and to increase the performance of muscle groups. There are two different types of EMS application. One is for targeted strengthening of muscles (activating application), and the other is to achieve a relaxing, restful effect (relaxing application through massage-like stimulation).

5. PACKAGE CONTENTS

1x EM 46
2x connection cables
8x adhesive electrodes 1.8 x 1.8 in (45 x 45 mm)
1x belt clip
3x AAA batteries
1x instructions for use


6. PARTS AND CONTROLS

Buttons:

1. ON/OFF button / **ENTER** button ⏻
2. Return button ⮌/ Favorite program ♥
3. Setting buttons **</>** (for selecting the category and program)
4. Intensity setting buttons **+ / –**
5. Key lock 🔒
6. Electrode connection, channels 1 and 2

Display (full screen):

7. Menu: (EMS) / (TENS acute pain) / (TENS chronic pain) / (RELAX / Massage-like sensation (MLS))
8. Favorite program set ♥
9. Key lock 🔒
10. Low battery level
11. Program number
12. Timer function (remaining time display)
13. Electrode positioning indicator
14. Pulse intensity, channel 1
15. Pulse intensity, channel 2


7. INITIAL SETUP

1. Remove the belt clip from the unit (if attached).
2. Press the battery compartment cover on the rear of the device and slide it downwards.
3. Insert three “AAA” alkaline batteries in accordance with the polarity markings inside the compartment.
4. Reattach the battery compartment cover. (Fig. 1)
5. Attach the connection cable to the electrodes. (Fig. 2)


ⓘ The electrodes feature clip closures to make attachment particularly easy.

6. Insert the plugs of the cables into the sockets on the top of the device. (Fig. 3)
7. Attach the belt clip, if desired.
8. Do not pull, twist or make sharp bends in the cables. (Fig. 4)


(Fig. 1)


(Fig. 2)

(Fig. 3)

(Fig. 4)

⚠ WARNING

Peel off the protective film slowly and carefully. Make sure that the adhesive electrode is not damaged, as damage or unevenness may cause skin irritation.

8. TO USE

Starting use

Step 1: Place the electrodes on your target area (for positioning suggestions, see page section on “Electrode positioning”) and connect them to the device.

Step 2: Press the ON/OFF button ⏻ to switch the device on.

Step 3: Press the **</>** buttons to navigate through the (EMS) / (TENS) / (RELAX/MLS) menus and press the ⏻ button to confirm your selection. If you have selected the TENS category, press the **</>** buttons to select your preferred subcategory or .

Step 4: Use the **</>** setting buttons to select the program number you want and press the ⏻ button to confirm your selection. The figure on the display indicates which parts of the body each program was developed for.

Step 5: Use the **+/–** setting buttons to select your preferred pulse intensity. If you want to control the channels individually, use the **< />** setting buttons to select the relevant channel. A line below the intensity indicator shows which channel has been selected. Change its pulse intensity using **+/–**. The indicator for pulse intensity in the display changes accordingly.

Notice

- After 10 seconds of inactivity, the device’s illuminated display switches off automatically. You can reactivate it by pressing any button (except the ON/OFF button ⏻).
- The device switches itself off automatically if it is not used for 3 minutes (automatic switch-off). When switched on again, the LCD screen shows the category selection (EMS) / (TENS) / (RELAX/MLS) – the most recently used category will be illuminated.
- If you press a button that you should not press, a short beep will sound.

ⓘ GENERAL INFORMATION

Return ⮌

If you wish to return to the previous selection menu, press the Return button ⮌.

Key lock

Locks the buttons during stimulation to avoid them being pressed accidentally.

1. To activate, hold down the **button lock** button at the side for approx. 2 seconds until the symbol displays.
2. To deactivate, hold down the **button lock** button again for approx. 2 seconds until the symbol disappears.

Pause

To pause the stimulation at any time, press the ON/OFF button ⏻. To continue the stimulation, set your preferred pulse intensity again.

9. PROGRAM OVERVIEW

The digital TENS/EMS unit features a total of 22 programs that are tailored to the different areas of the body:

- 4 TENS programs for treating acute pain
- 5 TENS programs for treating chronic pain
- 8 EMS programs for muscle stimulation
- 5 RELAX / Massage-like sensation (MLS) programs for muscle relaxation

In all of the programs you can set the pulse intensity of the two channels individually.

The figure on the display is intended as an initial aid to help you position the electrodes.

EMS Program Table

Prog. no.	Effective application areas, indications	Running time (min)	Possible electrode positioning
1	Lower back and back of thigh	40	09, 10, 23
2	Lower arms and calf	40	15, 16, 26, 27
3	Shoulder, elbow, wrist, knee and ankle	40	11, 12, 17, 18, 25, 28
4	Front of thigh	15	24
5	Bottom	15	22
6	Shoulders, upper back and upper arms	15	01–08, 11-14
7	Abdomen	15	19–21
8	Lower abdomen	20	21

TENS Acute Pain Program Table

Prog. no.	Effective application areas, indications	Running time (min)	Possible electrode positioning
1	Abdomen	20	19–20
2	Back of thigh	20	23
3	Shoulder, upper and lower back and calf	20	01–12, 26, 27
4	Arms, front of thigh, shoulder, elbow, wrist, knee and ankle	20	13-16, 23

TENS Chronic Pain Program Table

Prog. no.	Effective application areas, indications	Running time (min)	Possible electrode positioning
1	Shoulder, arms, elbow, wrist, knee and ankle	30	11–18, 25, 28
2	Calf	30	26, 27
3	Abdomen and front of thigh	30	19–20, 24
4	Shoulder and upper back	30	01–08, 11, 12
5	Lower back and back of thigh	30	09, 10, 23

RELAX/MLS Program Table

Prog. no.	Effective application areas, indications	Running time (min)	Possible electrode positioning
1	Abdomen and front of thigh	25	19–21, 24
2	Shoulder, upper and lower back and back of thigh	25	01–12, 23
3	Arms	25	13–16
4	Calf	25	26, 27
5	Shoulder, elbow, wrist, knee and ankle	25	11, 12, 17, 18, 25, 28

10. FAVORITE PROGRAM

With the favorite program, you can define a favorite from the 22 TENS / EMS / RELAX/MLS programs.
This makes it easier and quicker for you to access your favorite program.
If you have set a favorite program, when you switch on the device the favorite program will be automatically accessed. You can then start stimulation directly in your favorite program. A heart ♥ on the display indicates that you are in the favorite program.
The choice of favorite program may be based on your personal preference or, for example, the advice of your doctor.

Setting the favorite program

1. Select your desired program from the 22 available programs.
2. To set the selected program as a favorite, press and hold the favorite button ♥ for 2 seconds.
3. The saving of the favorite program is confirmed by a short acoustic signal. The ♥ symbol also appears on the display. This shows you that you are in the favorite program. When the device is switched on again, your favorite program is accessed directly.

 Note: Even if you have set a favorite program, you can access the program selection by pressing the Return button ⮌. You can replace your favorite program with another one at any time.

Deleting the favorite program

To delete the favorite program, press and hold the favorite button ♥ for approx. 2 seconds. The deletion of the favorite program is confirmed by a short acoustic signal. The previously displayed favorite program symbol ♥ goes out.

11. THERAPY MEMORY

Your device records the treatment time in the therapy memory. This allows you to record for how long you have been using the device, which may be helpful when consulting your doctor.


Accessing the therapy memory

To access the therapy memory, switch on the device using the ON/OFF button and, if you have activated a favorite program, press the Return button ⮌ to go back to the first level of the menu. Press and hold the **>** button for 3 seconds.
The elapsed treatment time appears in the display. The left two numbers stand for hours; the minutes are shown on the right. To exit the therapy memory display, press the Return button ⮌.

Resetting the therapy memory

To reset the treatment time memory (therapy memory) to 00, press and hold the **–** button for 3 seconds.

12. ELECTRODE POSITIONING

The figure on the display is intended as an initial aid to help you position the electrodes.

⚠ WARNING

Do not apply the electrodes to the pectoral muscles (chest). The device must never be used close to the heart. The stimulation electrodes must not be placed on any part of the front ribcage (where the ribs and breastbone are located), especially not on the two large pectorals. This can increase the risk of ventricular fibrillation and induce cardiac arrest.


It is critical to position the electrodes properly. Consult your doctor to establish the ideal electrode positions for your intended application area. The above chart is intended as an initial aid to help position the electrodes.

Electrode Distance

The larger the distance between electrodes, the larger the stimulated tissue volume. This applies to the area and the depth of the tissue volume. At the same time, however, the stimulation intensity of the tissue falls the further the electrodes are apart. As a result, greater distances between electrodes mean a larger tissue volume is stimulated, but less intensively. Consequently, you must increase the impulse intensity to boost stimulation.

- Proper distance between electrodes: 2-6 inches (5-15 cm)
- At distances less than two inches (5 cm), the unit primarily stimulates surface structures intensively.
- At distances more than 6 inches (15 cm), large areas and deep structures are stimulated very weakly.

Relationship between electrodes and muscle fiber structures

Adapt the current flow direction to the fiber structure of the muscle, according to the muscle layer you would like to treat. If you are targeting surface muscles, position the electrodes parallel to the fiber structure (A–B / C–D) and if you are targeting deeper layers of tissue, position the electrodes across the fiber structure. You can do this by positioning electrodes as crosses (i.e. diagonally), such as A–D/B–C.

NOTE: As part of pain relief treatment (TENS) using separately adjustable channels and two adhesive electrodes each, either position both electrodes of a channel so the area affected by pain is between the electrodes or position one electrode directly on the area affected by pain and the other electrode at least 1-1½ inches (2-3 cm) away.

You may use the electrodes of the second channel to simultaneously treat additional areas affected by pain or use them in conjunction with the electrodes of the first channel to restrict the area affected by pain (position electrodes opposite).
In this case, position electrodes as crosses.

Tips:


- Tip for the RELAX/MLS function: Always use at least 4 electrodes for optimum treatment.
- Use the electrodes on skin that is clean and preferably free from hair and grease in order to prolong the life of the electrodes. For best results, clean the skin with water and remove hair prior to treatment.


– If an electrode comes loose during use, the impulse intensity of the corresponding channel is reduced to the lowest level. Apply the electrode again and reset your desired impulse intensity.

13. ELECTRIC CURRENT PARAMETERS

Electrostimulation devices operate with the following electrical current parameters, which may affect the stimulation results differently depending on the setting:


Pulse shape


This describes the time function of the stimulating current.
A distinction is made between monophasic and biphasic pulse currents. In monophasic pulse currents the current flows in one direction, and in biphasic pulse currents the electrical pulse alternates its direction.
The device only provides biphasic pulse currents, as these relieve muscles, cause little muscle fatigue and are safe for use.

Pulse frequency


The frequency indicates the number of individual pulses per second and is displayed in Hz (Hertz). It can be calculated by determining the cyclic value for the time period. The relevant favorably determines which types of muscle fibres react most favourably. Slow-reaction fibres react more easily to lower pulse frequencies of up to 15 Hz, whereas rapid-reaction fibres only respond from approximately 35 Hz up.
Pulses of approx. 45-70 Hz are linked with constant tension in the muscles and quicker fatigue. Higher pulse frequencies are therefore better to use for speed strength and maximum power training.

Pulse width

This indicates the duration of an individual pulse in microseconds. Among other aspects, the pulse width therefore determines the penetration depth of the electricity, where usually larger muscle mass requires a larger pulse width.

Pulse intensity

Adjusting the intensity level depends on the individual sensitivity of each user and is determined by a variety of variables, such as site of application, blood supply to the skin, skin thickness and the quality of the electrode contact. The setting should be effective but must never cause any unpleasant sensation, such as pain at the site of application. While a gentle tingling indicates sufficient stimulation energy levels, any setting that causes pain should be avoided.
If using the device over an extended period, you may need to adjust the intensity level as your muscles start to adapt to the pulse intensity at the site of application.

Cycled pulse parameter variation

In many cases it is necessary to cover the overall tissue structure at the site of application by applying several pulse parameters. In this device, this is achieved by the provided programs, which automatically make a cyclical pulse parameter change. This also prevents individual muscle groups at the site of application being affected by fatigue.

14. CARE, MAINTENANCE AND DISPOSAL

Adhesive electrodes

- To ensure that the adhesive electrodes remain adhesive for as long as possible, clean them carefully with a damp, lint-free cloth and pat dry with a lint-free cloth. Before cleaning with water, remove the connection cables from the electrodes.
- Following treatment, stick the electrodes back onto the carrier foil.

Cleaning the device

- Clean the device after use with a soft, slightly damp cloth. If it is very dirty, you can also moisten the cloth with a mild soapy solution. Ensure that no water enters the device.
- Do not use any chemical or abrasive cleaning agents for cleaning.

Storage

- Do not make sharp bends in the connection cables or the electrodes.
- After use, stick the electrodes back onto the carrier foil for the gel electrodes.
- Store the device in a cool, well-ventilated location.
- Never place any heavy objects on the device.
- If the device is not going to be used for a long period of time, remove the batteries from the battery compartment.

Disposal

Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

15. TROUBLESHOOTING GUIDE

Problem: The device does not switch on.
Solutions:
(1) Ensure batteries have been inserted correctly and are in contact with the terminals.
(2) Replace the batteries if necessary.
(3) Contact Customer Services.

Problem: Electrodes do not adhere to the body.
Solutions:
(1) Clean the adhesive surface of the electrodes using a damp, lint-free cloth. If they still do not stick securely, the electrodes must be replaced.
(2) Clean the skin prior to any application; do not use skincare lotions or oils prior to treatment. Shaving may increase the life of electrodes.

Problem: There is no noticeable stimulation.
Solutions:
(1) Press the ON/OFF button ⏻ to interrupt the program. Check that the connection cables are correctly connected to the electrodes. Ensure that both electrodes of the channel you are using are in firm contact with the treatment area. Check the electrode positions and ensure that the adhesive electrodes do not overlap.
(2) Ensure the connection plug is firmly connected to the device.
(3) Gradually increase the intensity of the affected channel.
(4) If you still cannot feel any stimulation, replace the batteries.

Problem: The low battery indicator displays.
Solution:
Replace all batteries.

Problem: Unpleasant sensation near the electrodes.
Solutions:
(1) Reduce the intensity of the affected channel.
(2) The electrodes are not positioned correctly. Check and reposition.
(3) The electrodes are worn. This may cause irritated skin because of uneven current across the area. Replace the electrodes.

Problem: Skin in the treatment area turns red.
Solutions:
Immediately stop treatment and wait for the skin color to return to normal. If the redness is directly under the electrode and goes away after a brief time, this is normal and is caused by increased blood flow. Consult your doctor before continuing if the skin irritation persists and is accompanied by an itching sensation or inflammation, as this may be an allergic reaction to the adhesive.

Problem: A recurring beep sounds, and the electrode intensity cannot be increased.
Solution:
(1) The adhesive electrodes are not correctly applied to the skin. Check the adhesion and re-attach them if necessary.
(2) The adhesive electrodes may be worn out and are no longer able to conduct the current. Please replace them.

16. TECHNICAL SPECIFICATIONS

Model	EM 46
Output waveform	Biphasic rectangular and needle pulses
Pulse length	50–400 µs
Pulse frequency	2–120 Hz
Output voltage	max. 100 Vpp (at 500 Ohm)
Output current	max. 200 mApp (at 500 Ohm)
Power supply	3 x AAA alkaline batteries (LR03)
Treatment time	15–40 minutes
Intensity levels	adjustable from 0 to 50
Electrodes used	Self-adhesive gel electrodes 1.8 x 1.8 in (45 x 45 mm)
Operating conditions	41°F to104°F (5°C to 40°C) at a relative air humidity of 15–90%
Storage conditions	32°F to 104°F (0°C to 40°C) at a relative air humidity of ≤90%
Transport conditions	-13°F to 158°F (-25°C to 70°C) at a relative air humidity of ≤90%
Dimensions	Approx. 5.2 × 2.6 × 1 in (132 x 66 x 26 mm)
Weight	Approx. 3.67 oz (104 g)
Altitude limit for use	3,000 m
Maximum tolerated atmospheric pressure	700 hPa to 1060 hPa

17. REPLACEMENT AND CONSUMABLE PART

You can obtain the following replacement parts directly from www.shop-beurer.com:

Designation	Item number
20 x adhesive electrodes 1.8 x 1.8 in (45 x 45 mm)	110.327
Connection cables	110.392

18. FCC COMPLIANCE INFORMATION

Digital TENS/EMS device EM 46

Responsible Party – U.S. Contact Information
Beurer North America LP, 1 Oakwood Boulevard, Suite 255,
Hollywood, FL 33020, United States
1-800-536-0366
info@beurer.com

FCC Compliance Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

MODIFICATION: Any changes or modifications not expressly approved by the grantee of this device could void the user's authority to operate the device.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

19. WARRANTY

5-Year Limited Warranty for Original Purchaser

Your Beurer Digital TENS/EMS device, Model EM 46, is warranted to be free from defects in materials and workmanship for a period of five (5) years from the date of purchase under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the Beurer Digital TENS/EMS device, Model EM 46, without additional charge, for any part or parts covered by this written warranty. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.

For warranty service, contact our customer service department at 1-800-536-0366 or info@beurer.com with a description of the problem. If the issue is deemed to fall within the scope of the 5-Year Limited Warranty, you will be asked to mail the product (at your cost) in its original package along with proof of purchase, your name, address, and phone number. If the issue is not deemed eligible under warranty, we will provide a quotation for repair or replacement, including return shipping fees.

Third-Party Repairs

This warranty complies with the federal Magnuson–Moss Warranty Act. You may use third-party parts or services for repairs without voiding this warranty, unless the repair itself causes damage to the

product. We encourage you to retain your proof of purchase and contact our customer service team with any questions.

Unauthorized Sellers

This warranty applies only to products purchased from Beurer or authorized retailers and distributors. Products purchased from unauthorized sellers — including third-party platforms or online marketplaces (such as eBay or unverified Amazon sellers) — are not covered by this warranty. This policy is in place to protect customers from counterfeit, altered, or used products that may not meet Beurer's quality standards.

If you're unsure whether your seller is authorized, please contact our support team.

Exclusions

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessories; alterations; improper installation; misapplication; lack of reasonable care; unauthorized repairs or modifications; improper electrical/power supply; worn batteries; normal wear and tear; power loss; dropped product; failure to follow instructions; lack of recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions.

This warranty is void if the product is used in a commercial or business environment.

Beurer's maximum liability under this warranty is limited to the actual purchase price paid, as confirmed by proof of purchase. The warranty is effective only in the country where the product was purchased and used.

Legal Notice

This warranty gives you specific legal rights, and you may have other rights which vary from state to state.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE LENGTH OF THIS WRITTEN WARRANTY.

IN NO EVENT SHALL BEURER BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING FROM USE OR INABILITY TO USE THE PRODUCT.

Beurer does not authorize any third party — including retailers or resellers — to modify the terms of this warranty.

This warranty does not apply to opened, used, repackaged, or resold products, including those sold through online auctions or unauthorized channels.

For more information regarding our product line in the USA, please visit: www.shop-beurer.com

For Canada only:
Manufacturer:
Beurer GmbH
Söflinger Straße 218
89077 Ulm, Germany

Distributed by:
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020, USA
www.shop-beurer.com

Questions or comments? Call our US-based customer service toll free at 1-800-536-0366.

Assembled in Indonesia

Subject to errors and changes

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The device may not be fully usable in the presence of electromagnetic interference. This could result in issues such as error messages or the failure of the display/device.
- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner outlined above, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of components other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device‘s resistance to electromagnetic interference; this can result in faulty operation.
- Keep portable RF communication devices (including peripheral equipment such as antenna cables or external antennas) at least 2.54 in (30 cm) away from all device parts, including all cables included in delivery. Failure to comply with the above could impair the performance of the device.
- The use of components other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device’s resistance to electromagnetic interference; this can result in faulty operation.
- Keep portable RF communication devices (including peripheral equipment such as antenna cables or external antennas) at least 30 cm (12 inches) away from all device parts, including all cables included in delivery. Failure to comply with the above could impair the performance of the device.

Table 1

Guidance and manufacturer’s declaration - electromagnetic emissions

Emissions test	Compliance
RF emissions, CISPR 11	Group 1,
RF emissions, CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/ flicker emissions, IEC 61000-3-3	Not applicable

Table 2

Guidance and manufacturer’s declaration - electromagnetic Immunity

Immunity Test	IEC 60601 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±4 kV ±8 kV contact ±2 kV ±4 kV ±8 kV ±15 kV air	±4 kV ±8 kV contact ±2 kV ±4 kV ±8 kV ±15 kV air
Electrical fast transient/ burst IEC 61000-4-4	Not applicable	Not applicable
Surge IEC 61000-4-5	Not applicable	Not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Not applicable	Not applicable

Immunity Test	IEC 60601 Test level	Compliance level
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3V(0.15-80 MHz); 6V in ISM and amateur radio bands between 0.15-80 MHz	3V(0.15-80 MHz); 6V in ISM and amateur radio bands between 0.15-80 MHz
Radiated RF IEC61000-4-3	10V/m 80-2700 MHz	10V/m

Table 3

Manufacturer's declaration - electromagnetic Immunity to proximity fields from RF wireless communications equipment

	IEC 60601 Test level				Compliance Level (V/m)
	Test frequency (MHz)	**Modulation**	**Maximum Power (W)**	**IMMUNITY TEST LEVEL (V/m)**	
Radiated RF IEC61000-4-3	385	Pulse modulation 18 Hz	1,8	27	27
	450	FM c) ± 5 kHz deviation 1 kHz sine	2	28	28
	710	Pulse modulation 217 Hz	0.2	9	9
	745				
	780				
	810	Pulse modulation 18 Hz	2	28	28
	870				
	930				
	1720	Pulse modulation 217 Hz	2	28	28
	1845				
	1970				
	2450	Pulse modulation 217 Hz	2	28	28
	5240	Pulse modulation 217 Hz	0.2	9	9
	5500				
	5785				

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