MG 35 Deep Roll

ENGLISH

Deep Roll

Instruction for Use.....3

ESPAÑOL

Rodillo de masaje Deep Roll

Instrucciones de uso.....26

FRANÇAIS

Appareil de massage profondeur

Mode d'emploi53

Distributed by/
Distribuido por/
Distribué par :
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020, USA
www.beurer.com

Questions or Comments? Call our US-based customer service toll free at 1-800-536-0366.

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READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT

Keep this manual in a safe location for future reference

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ENGLISH

Contents

1. Important safety instructions 5
2. Introduction 12
3. Package Contents 12
4. Parts and Controls 13
5. To Use 14
6. Care, Maintenance and Disposal 19
7. Technical Specifications 19
8. FCC Compliance Information 20
9. Warranty 21

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE
If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact customer service. See warranty for service contact.

Dear Valued Customer,
Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for application in the areas of heating, weight management, blood pressure, body temperature, pulse, gentle therapy, massage, and beauty.

Register your product today to receive further benefits:
- added protection for your purchase
- health and wellness tips
- new product information
- continual outstanding customer service and support.

Visit www.registerbeurer.com or contact customer service at 1-800-536-0366. Please read these instructions for use carefully and keep them for later use. Make them accessible to other users and observe the information they contain.

With kind regards, Your Beurer team


1. Important safety instructions

Signs and symbols

Whenever used, the following signs identify safety and property damage messages and designate a level of hazard or seriousness.

READ THIS ENTIRE MANUAL, THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

⚠	This is the safety alert symbol. It alerts you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
WARNING	Indicates a potentially hazardous situation, which, if not avoided, could result in death or serious injury.

CAUTION	Indicates a potentially hazardous situation, which, if not avoided, may result in minor/ moderate injury.
NOTICE	Addresses practices not related to personal injury, such as product and/or property damage.
 PIN1 “+” PIN4 “-”	Polarity

WARNING – To reduce the risk of fire, electric shock, or serious personal injury:

- The device may only be used for the purpose for which it is designed and in the manner specified in the instructions for use. Any improper use may be dangerous.
- The manufacturer is not liable for damage resulting from improper or incorrect use.
- In the event of a malfunction, shut the massager off immediately and disconnect it from the powersupply.

- Do not insert any objects into the opening in the device or into the moving parts. Ensure that the moving parts can move freely at all times.
- The device is suitable only for massage treatment of the thighs, calves, shoulders, arms and back region (excluding the spine). The device is not a medical device; rather it is a massage treatment device.
- Do not use the massager on the head. Keep hair, jewelry, and loose clothing away from the massage head.
- Do not use the massager on the spine
- The device is intended only for domestic/private use, not for medical or commercial use.
- Only use the device with the supplied accessories.
- Keep packaging material away from children. There is the risk of suffocation.
- Never use underneath a cover, such as a blanket, pillow, etc.
- Never use near gasoline or other highly flammable substances.
- Only use indoors but never in the bath or sauna.
- Keep device away from water or other liquids.
- Never reach for a device that is submerged in water.
- Never leave this device unattended, especially if children are nearby

- Do not use massager on children, animals, or on persons suffering from a herniated disk or other spinal problems.
- Children must not play with the device.
- Cleaning and user maintenance must not be performed by children unless supervised.
- Under no circumstances should you open or repair the device yourself, as faultless functionality can no longer be guaranteed thereafter. Failure to comply will invalidate the warranty.
- Only specialist personnel may perform repairs on electrical devices. Improper repairs may subject users to considerable danger. For repairs, please contact Customer Services or an authorised retailer.
- Do not use the massager if one or several of the following warnings applies to you. If you are unsure whether the device is suitable for you, please consult your doctor.
 - On animals
 - If you suffer from a medical abnormality or injury in the area to be massaged (e.g. herniated disc, open wounds)
 - Under any circumstances for massages in the area around the heart if you wear a pacemaker. If this is the case, consult your doctor regarding other parts of the body first

 - On the area around the heart
 - On swollen, burnt, inflamed or damaged skin and body parts of any kind
 - On cuts, capillaries, varicose veins, acne, couperose herpes, or other skin diseases
 - During pregnancy
 - On the face (eyes), on the larynx or other particularly sensitive body parts
 - While sleeping
 - After the consumption of medicines or alcohol (reduced cognition!)
 - In a vehicle
 - In the pubic area
- Consult your doctor before using the massager, especially
 - If you suffer from a severe illness or are recovering from an operation
 - If you have a pacemaker, implant or other aids
 - If you suffer from thrombosis
 - If you suffer from diabetes
 - If you are from suffering pain, the cause of which is unknown.

– Any of the following: Cold accompanied by a fever, varicose veins, phlebitis, jaundice, nervous system diseases (e.g. sciatica) or acute inflammation.

CAUTION:

- Shut the device off after every use and before cleaning.
- Before use, ensure that the device and accessories show no visible damage. If you have any doubts, do not use the device and contact your retailer or the specified Customer Service address.
- Never open the housing.
- Keep the device away from pointed or sharp objects.
- If the device has been dropped or has suffered any other damage, it must no longer be used

Notes on electromagnetic compatibility

- The device is suitable for use in all environments listed in these instructions for use, including domestic environments.
- The use of the device may be limited in the presence of electromagnetic disturbances. This could result in issues such as error messages or the failure of the display/device.

- Avoid using this device directly next to other devices or stacked on top of other devices, as this could lead to faulty operation. If, however, it is necessary to use the device in the manner stated, this device as well as the other devices must be monitored to ensure they are working properly.
- The use of accessories other than those specified or provided by the manufacturer of this device can lead to an increase in electromagnetic emissions or a decrease in the device's electromagnetic immunity; this can result in faulty operation.
- Failure to comply with the above can Impair the performance of the device.

Battery Handling Safety Precautions

- Protect batteries from excessive heat.
- Do not dispose of batteries in fire, batteries may explode or leak.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Keep batteries away from children and pets. Batteries may be harmful if swallowed. Should a child or pet swallow a battery, seek medical assistance immediately.
- Do not disassemble, split or crush the batteries.
- Only use chargers specified in the instructions for use.

- Only charge the device using the supplied AC adapter.
- In order to achieve as long a battery service life as possible, fully charge the battery at least 2 times a year.
- Always switch the device off in advance of charging.
- Batteries must be charged correctly prior to use. The instructions from the manufacturer and the specifications in these instructions for use regarding correct charging must be observed at all times.

2. Introduction

This Deep Roll enables you to target tense muscle groups for a precise massage. It can also help loosen fascial tissue and break down fasciae agglutinations. The Deep Roll is comfortable to use thanks to its soft-touch surface. It has three independently moving sections to adapt to your body movements. It also features three different intensity levels and a preset massage sequence.

3. Package Contents

1 x Deep Roll MG 35
1 x USB charging cable
1 x Instruction Manual

4. Parts and Controls


1 ON/OFF button	
2 LED display (3 LEDs)	
Left LED lights up	= Low intensity level
Centre LED lights up	= Medium intensity level
Right LED lights up	= High intensity level
All LEDs light up	= Pre-set massage mode
3 Connection for USB charging cable	

5. To Use


Charging the Battery

Before using the Deep Roll for the first time, charge it first for a minimum of 2 ½ hours. Proceed as follows:

1. Connect the supplied micro USB cable to the USB port on the side of the Deep Roll opposite the Power Button.

2. Attach the larger USB end of the cable to a computer or other suitable 5VDC charging source.

3. Attach the larger USB end of the cable to a computer or other suitable 5VDC charging source.


CAUTION: The massage should always feel pleasant and relaxing. If you feel pain or discomfort, stop the massage, change your position, or decrease the contact pressure. A slight hematoma may well occur – this is a positive effect as the use of the Deep Roll can selectively encourage your blood flow. In isolated cases, areas of redness may also occur on the skin. If you experience excessive skin irritation, stop using the device and consult a doctor.

NOTE: We recommend massaging each area (lower body, mid-body and upper body) three times a week for 15 minutes each time. When carrying out a session, we recommend massaging all three areas. However, you may also choose to focus on one of the above areas. Carry out the massage consciously and slowly.

Switch on the device by pressing the Power Button once; the Deep Roll starts to vibrate at the low intensity level. Press the Power Button again to increase intensity as desired. Pressing the Power Button four times will cause the Deep Roll to enter a preset massage sequence in which it randomly starts and stops in various intensity levels and durations. To shut the device off at any time, press and hold the Power Button for about two seconds.


Treatment Suggestions

Calf


Sit on the floor and place one leg on the Deep Roll so that it is positioned on the bottom half of the calf. Support yourself with your hands while leaning back and lift your hips. Tense your stomach and bottom muscles and shift your lower body forward and backward.

Rear Thigh/Glutes


Lie on your side on the floor and stretch both legs straight out. Bend one arm and support yourself on the floor with your elbow while keeping your other arm relaxed by your side. Place the Deep Roll under your hips. Now roll over the Deep Roll with slow movements. Vary the exercise by supporting yourself with both arms while leaning back; then roll your bottom over the Deep Roll with gentle pressing movements.

Thigh


Sit upright on the floor or on a chair and hold the Deep Roll by the two outer sections. Roll it over your thigh slowly, applying pressure as desired. The middle section of the Deep Roll will treat the thigh muscle.

Arm/Shoulders

Lean one side against a wall while standing. Place the Deep Roll between your upper arm and the wall. Apply gentle pressure and move up and down with your legs so that the Deep Roll moves slowly along your arm. Increase the distance between your legs and the wall to increase pressure as desired. Next, slowly turn your back to the wall so that the Deep Roll moves along your shoulder area.

Back

Lie on your back and place the Deep Roll to the right or left of the spine between your back and the floor. Bend your legs and move slowly back and forth to roll the Deep Roll along your back.

CAUTION: Do not roll the Deep Roll directly over bones.

6. Care, Maintenance and Disposal

NOTICE:

- Never let water enter the device.
- Clean the device with a slightly damp cloth. Do not use abrasive cleaning products or stiff brushes.
- Make sure no liquids enter the device.
- Ensure device is fully dry before use.
- Do not clean the device in a dishwasher.

Disposal

Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

7. Technical Specifications

Type	MG 35
Dimensions	6.1 x 3.4 x 3.4 inches (15.6 x 8.7 x 8.7 cm)
Weight	26.3 ounces (746 g)
Internal Battery Pack	1800 mAh 3.7 V Lithium-ion

8. FCC Compliance Information

Deep Roll MG 35
Responsible Party – U.S. Contact Information
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020, United States
1-800-536-0366
info@beurer.com

FCC Compliance Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
Changes or modifications to the product not expressly approved by the manufacturer could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide

reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

9. Warranty

Limited Lifetime Warranty For Original Purchaser

Your Beurer Deep Roll, model MG 35, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service.

This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.
We will, at our option, repair or replace the Beurer Deep Roll, model MG 35, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.
For warranty service contact our customer service department at 1-800-536-0366 or at info@beurer.com to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.
This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessory; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized repairs or modifications; improper use of electrical/power supply;

old worn batteries; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Beurer. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of Beurer under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL BEURER BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

Beurer does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate Beurer in any way beyond the terms set forth herein.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or

products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of Beurer.
This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.
For more information regarding our product line in the USA, please visit: www.beurer.com

Distributed by:
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020, USA
www.beurer.com

Questions or comments? Call our US-based customer service toll free at 1-800-536-0366.

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