beurer

TB 10

ENGLISH Electric Toothbrush Set for Kids
Instructions for use 3

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READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT

Keep this manual in a safe location for future reference

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ENGLISH

Contents

1. Important Safety Instructions 4
2. Package Contents 6
3. Parts and Controls 7
4. Introduction 7
5. To Use 8
6. Care, Maintenance and Disposal 9
7. Technical Specifications 10
8. Warranty 10

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact Beurer customer service. See warranty for service contact.

Dear Valued Customer,
Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for application in the areas of heating, weight management, blood pressure, body temperature, pulse, gentle therapy, massage, and beauty.

Register your product today to receive further benefits:
- added protection for your purchase
- health and wellness tips
- new product information
- continual outstanding customer service and support.

Visit www.registerbeurer.com or contact customer service at **1-800-536-0366.**

Please read these instructions for use carefully and keep them for later use. Make them accessible to other users and observe the information they contain.

With kind regards,
Your Beurer Team

1. IMPORTANT SAFETY INSTRUCTIONS

Signs and symbols

Whenever used, the following signs identify safety and property damage messages and designate a level of hazard or seriousness.

READ ALL INSTRUCTIONS BEFORE USING

⚠	This is the safety alert symbol. It alerts you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
WARNING	Indicates a potentially hazardous situation, which, if not avoided, could result in death or serious injury.
CAUTION	Indicates a potentially hazardous situation, which, if not avoided, may result in minor/moderate injury.
NOTICE	Addresses practices not related to personal injury, such as product and/or property damage.

⚠ WARNING – To reduce the risk of fire, electric shock, or serious personal injury:

- This product is not a toy; adult supervision is needed whenever children use this product to help prevent a choking hazard or misuse. Do not let children play with this product.
- Cleaning and user maintenance must not be performed by children unless supervised.
- Usage by children under age 3 is not recommended.

- Do not disassemble, alter, or attempt to repair this device.
- Do not use this device while bathing or showering.
- Do not immerse device (handle) in water or other liquids.
- This product is not intended for use by persons with reduced physical, sensory, or mental capabilities, or lack of experience and knowledge, unless they are supervised or have been given instruction by a person responsible for their safety.
- If excessive discomfort or pain is experienced during use, discontinue use and contact your dentist.
- Shut the product off immediately if it malfunctions during use, before maintenance, when changing attachments, and when finished using.

⚠ CAUTION – To reduce the risk of personal injury or product/property damage:

- This device is intended solely for household use in oral cleaning. Do not use this device for any other purpose.
- Discontinue use if any parts of this product appear damaged. This product has no user serviceable parts; if it becomes inoperable or needs repair, contact Customer Service for warranty coverage.
- A brush head is intended for a single user; do not allow multiple users to use one brush head.
- Do not use attachments other than the ones included with this product or recommended in this manual.
- Do not press brush head against teeth and gums too firmly.

Battery Handling Safety Precautions

- Use only the size and type of batteries specified.
- Be sure to follow the correct polarity when installing the batteries. Reversed batteries may cause damage to the device.
- Do not mix different types of batteries together (e.g. Alkaline and Carbon-zinc or rechargeable batteries) or old batteries with fresh ones. Always replace batteries as a simultaneous set.
- If the batteries in the device are depleted or the device will not be used for a long period of time, remove the batteries to prevent damage or injury from possible battery leakage.


- Do not try to recharge batteries not intended to be recharged; they can overheat and rupture (follow battery manufacturer's directions.)
- Do not dispose of batteries in fire, batteries may explode or leak.
- Clean the battery contacts and also those of the device prior to battery installation.
- Remove discharged batteries from the product and dispose/recycle in compliance with all applicable laws.
- Keep batteries away from children and pets. Batteries may be harmful if swallowed. Should a child or pet swallow a battery, seek medical assistance immediately.

2. PACKAGE CONTENTS

Check that the exterior of the delivery cardboard packaging is intact and make sure that all components are present. Before use, ensure that the device and accessories show no visible damage and all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Service address.

1 x Electric Toothbrush Handle
2 x Brush Heads
1 x Animal Protective Cap
1 x Rinse Cup
1 x Sand Timer
1 x Storage Base
1 x Brushing Chart
1 x Instruction Manual

3. PARTS AND CONTROLS

1	Animal Protective Cap	2	Brush Head
3	Electric Toothbrush Handle	3a	ON/OFF Button
3b	Battery Compartment	4	Rinse Cup
5	Storage Base	6	Replacement Brush Head
7	Sand Timer		

4. INTRODUCTION

Only use the device for cleaning teeth on humans. Do not use the device on animals! The device is not suitable for persons under 3 years of age! Children should be supervised during use. The device is intended only for the purpose described in these instructions for use. The manufacturer cannot be held liable for damage caused by improper or careless use.

5. TO USE


Inserting/Replacing Batteries

1. Open the battery compartment by pulling the cover and insert two 1.5 V AA size alkaline batteries into the battery compartment according to the polarity markings.
2. Push the battery compartment cover with batteries back onto the handle.
3. Replace the batteries as soon as you notice brush head starts to turn slowly.

Attaching/Removing a Brush Head

NOTE: This product includes two brush heads. Both are identical except for the color for better identification.

1. After selecting a brush head, press it down onto the handle.
2. To remove the brush head, pull it away from the handle.


We recommend replacing brush heads after three months of usage or when brush head bristles become worn or flatten out.

Brushing

1. Wet bristles and apply toothpaste.
2. Turn over the 2 minutes sand timer.
3. To avoid splashing, hold the bristles against teeth at a slight angle and before pressing the ON/OFF button to activate the brush head spinning motion. Make sure to keep the mouth closed during brushing to prevent toothpaste splatter.

4. While brushing, move the brush head along the teeth and gums while applying light pressure. Make sure your child reaches all 4 jaw quadrants.
5. After the timer runs out, press the On/Off button to stop the spinning motion.
6. Rinse your mouth with water from the rinse cup.

NOTE: Gums may become sensitive after the first few uses if they are not accustomed to an electric toothbrush. This is normal and should subside with further usage. If sensitivity or discomfort persist, contact your dentist.

Brushing Chart

A Weekly Brushing Chart is included with this product to help keep track of regular brushing sessions. Additional charts can be downloaded from www.shop-beurer.com.

6. CARE, MAINTENANCE AND DISPOSAL

Cleaning

- Remove the brush head from handle and rinse it thoroughly after each use under clean running water.
- Clean the connection between the brush head and the handle.
- The device should be cleaned from time to time. Do not use any abrasive cleaning products and never place any parts in a dishwasher.
- Clean handle periodically with a moist cloth and a mild soap solution.
- Press the brush head back onto the handle.
- Put the handle into the base for storage and place the animal protective cover onto the brush head when it is dry.

- When not using this product for an extended period, remove the batteries.

Disposal

- Remove the batteries before disposing of the product.

Observe the local regulations for material disposal. Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

7. TECHNICAL SPECIFICATIONS

Type	TB 10
Size	Monkey: 235.9 x 160 x 60 mm Panda: 235.5 x 160 x 60 mm Tiger: 235.2 x 160 x 60 mm
Weight	0.37 lb (170 g)
Operating Time	120 Minutes
Operating Temperature Range	14° - 104° F (-10° - 40° C)
Storage Temperature Range	14° - 131° F (-10° - 55° C)

8. WARRANTY

Limited Lifetime Warranty For Original Purchaser
Your Beurer Electric Toothbrush Set for Kids, model TB 10, is warranted to be free from defects in materials and workmanship for the life of the product under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners.

We will, at our option, repair or replace the Beurer Electric Toothbrush Set for Kids, model TB 10, without additional charge, for any part or parts covered by these written warranties. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective

materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.

For warranty service contact our customer service department at 1-800-536-0366 or at info@beurer.com to provide a description of the problem. If the problem is deemed to be within the scope of the limited lifetime warranty, you will be asked to mail the product at your costs in its original package with proof of purchase, your name, address and phone number. If the problem is not deemed to be within the scope of the limited lifetime warranty, we will provide a quotation for repair respectively replacement and return shipping fee.

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessory; alteration to the product; improper installation; misapplication; lack of reasonable care with respect to the product; unauthorized repairs or modifications; improper use of electrical/power supply; old worn batteries; normal wear; loss of power; dropped product; malfunction or damage of an operating part as a result of failure to comply with instructions for use or to provide manufacturer's recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of Beurer. This warranty is void if the product is ever used in a commercial or business environment. The maximum liability of Beurer under this warranty is limited to the purchase price actually paid by the customer for the product covered by the warranty, as confirmed by proof of purchase, regardless of the amount of any other direct or indirect damage suffered by the customer.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adaptation to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE DURATION OF THIS APPLICABLE WRITTEN WARRANTY.
Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL BEURER BE LIABLE FOR ANY SPECIAL, INCIDENTAL, INDIRECT OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS, IMPLIED OR ANY OTHER THEORY OF LIABILITY, WHATSOEVER.
Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

Beurer does not authorize anyone, including, but not limited to, retailers, the subsequent consumer purchaser of the product from a retailer or remote purchaser, to obligate Beurer in any way beyond the terms set forth herein.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to sale of such products on Internet auction sites and/or products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate in connection with any products or parts thereof which are repaired, replaced, altered, or modified, without the prior explicitly written consent of Beurer.

This warranty gives you specific legal rights, and you may also have other rights which may vary from state to state.
For more information regarding our product line in the USA, please visit: www.shop-beurer.com

Distributed by:
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020, USA
www.shop-beurer.com

Questions or comments? Call our US-based customer service toll free at 1-800-536-0366.

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