BY 80

READ THIS MANUAL COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT
Keep this manual in a safe location for future reference

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ENGLISH **Baby scale**
Instruction manual 2

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Questions or comments? Call our US-based customer service toll free at 1-800-536-0366 or contact info@beurer.com

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ENGLISH

Contents

1. Signs and symbols 2
2. Important safety instructions 2
3. Package Contents 3
4. Battery Handling Safety Precautions 3
5. Before First Use 4
6. Parts and Controls 4
7. To Use 4
8. Care, Maintenance and Disposal 5
9. Technical Specifications 5
10. FCC Compliance Information 5
11. Warranty 6

SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

If the instruction manual is damaged or if you no longer have the instruction manual in your possession, please contact Beurer customer service. Please see warranty for service contact.

1. Signs and symbols

Whenever used, the following signs identify safety and property damage messages and designate a level of hazard or seriousness.

READ AND UNDERSTAND THIS ENTIRE MANUAL INCLUDING THE SAFETY SECTION AND ALL INSTRUCTIONS AND WARNINGS COMPLETELY AND CAREFULLY BEFORE USING THIS PRODUCT. FOLLOW ALL SAFETY INSTRUCTIONS AND WARNINGS TO AVOID HAZARDOUS SITUATIONS AND TO MAKE CORRECT USE OF THIS PRODUCT.

⚠	This is the safety alert symbol. It alerts you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.
⚠ **WARNING**	Indicates a hazardous situation, which, if not avoided, could result in death or serious injury.
⚠ **CAUTION**	Indicates a hazardous situation, which, if not avoided, may result in minor/moderate injury.
NOTICE	Addresses practices not related to personal injury, such as product and/or property damage.

2. Important safety instructions

⚠ WARNING - Safety notes

- Keep packaging material away from children (risk of suffocation).
- Swallowing batteries can be extremely dangerous. Keep the batteries and scale out of the reach of small children. Should a battery be swallowed, seek medical assistance immediately.
- Keep scale away from moisture, dust, chemicals, toiletries, liquid cosmetics, extreme hot and cold temperatures, and sources of heat (e.g., open flame, radiators).
- Prevent impact with hard objects.
- Do not press the button violently or with pointed objects.

⚠ IMPORTANT

- Do not place any object on the scale when it is not in use.
- Strong electromagnetic fields (e.g. cell phones) may impair the accuracy of the scale.
- Remove any transport locks.
- Never use harsh solvents or cleaning products!
- Do not submerge the baby scale in water.
- Do not rinse the baby scale under running water.

- Do not clean the baby scale in the dish washer.

General notes

- The scale has a capacity of 44 lb (20 kg). Measurements are displayed in 0.01 lb (5 g) divisions.

CAUTION

- Cleaning: You can clean the scale with a damp cloth and a mild soap if required.
- Never immerse the scale in water or rinse it under running water.
- The scale is not intended to be used for commercial purposes.
- Repairs may only be performed by an authorized service center. Contact Customer Service with any questions regarding operation, usage, or repair of this product.

Battery Handling Safety Precautions

WARNING - To reduce the risk of serious personal injury:

- If your skin or eyes come into contact with battery fluid, rinse the affected area with water and seek medical assistance.
- ⚠ **Choking hazard!** Small children may swallow and choke on batteries. Store batteries out of the reach of small children.
- Keep batteries away from children. Only adults should handle batteries.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect batteries from excessive heat.
- ⚠ **Risk of explosion!** Never throw batteries into a fire. Dispose or recycle spent/discharged batteries in compliance with all applicable laws.
- Follow the battery manufacturer's safety and usage instructions.
- Do not charge or short-circuit batteries.

CAUTION

- Always use the size and type of battery indicated.
- Do not mix a used battery with a new battery, or mix alkaline, standard (carbon-zinc), or rechargeable (nickel-cadmium) batteries.
- Insert the batteries observing the proper polarity (+/-) as indicated.
- Always replace all batteries at the same time.
- Do not use rechargeable batteries.
- Do not disassemble, open or crush batteries.

NOTICE

- Clean the battery contacts and those of the device prior to installation.
- Remove spent or discharged batteries from the product.
- Remove the batteries when not using the product for an extended period.

3. Package Contents

1 x Baby Scale
2 x 1.5 V AA alkaline Batteries
1 x Instruction for use

4. Battery Handling Safety Precautions

- Use only the size and type of batteries specified.
- Be sure to follow the correct polarity when installing the batteries. Reversed batteries may cause damage to the device.
- Do not mix different types of batteries together (e.g. Alkaline and Carbon-zinc carbon or rechargeable batteries) or old batteries with fresh ones. Always replace batteries as a simultaneous set.
- If the batteries in the device are depleted or the device will not be used for a long period of time, remove the batteries to prevent damage or injury from possible battery leakage.

- Do not try to recharge batteries not intended to be recharged; they can overheat and rupture (follow battery manufacturer's directions.)
- Do not dispose of batteries in fire, batteries may explode or leak.
- Clean the battery contacts and also those of the device prior to battery installation.
- Remove discharged batteries from the product and dispose/recycle in compliance with all applicable laws.
- Keep batteries away from children and pets. Batteries may be harmful if swallowed. Should a child or pet swallow a battery, seek medical assistance immediately.

5. Before First Use

Remove the Battery Compartment Cover located underneath the scale and insert two "AA" (LR06) size alkaline batteries (included) according to the polarity markings inside the compartment. Replace Battery Compartment Cover. When "Lo" appears in the LCD Display, the batteries need to be replaced.

6. Parts and Controls

1. Hold/Lock button
2. Display
3. Power/Tare button

7. To Use

WARNING

- **Never leave a child unattended while on the scale, even for a moment.**

1. Slide the Weight Selector Switch, located underneath the scale next to the Battery Compartment, to either "lb:oz" (pounds ounces), "lb" (pounds), or "kg" (kilograms).

NOTE: The default is set to "lb:oz" pounds ounces.
The scale has a capacity of 44 lbs (20 kg). Measurements are displayed in 0.01 lb (5 g) increments.

2. Place the scale on a flat, level, stable surface within your reach to ensure accurate measurement.
3. Press the Power Button to activate the scale. Weights are displayed in 0.01 lb (5g) increments.
4. When "0.000" appears on the display, the scale is ready for use. Carefully place the baby onto the scale, taking extra care with the baby's head.
5. This scale has both automatic and manual "HOLD" functions. The auto hold function displays the weight when it has stabilized enough to detect an exact value and remains displayed for 20 seconds.

 NOTE: If the baby is moving around too much, the auto hold will not function. In this case, press the HOLD Button when the baby is as calm as he/she will become. Within three seconds, the word "HOLD" stops flashing and the weight is displayed.
 If you do not want to place the baby directly on the weighing surface, but prefer to use a towel or blanket, place the towel/blanket on the scale first, and then press the Tare Button to "zero out" the display; "----" appears in the LCD Display and the weight resets to zero.

6. The LCD Display will shut off after about 20 seconds.

8. Care, Maintenance and Disposal

Care and Maintenance

- Clean the scale with a damp cloth and a little washing up liquid, if required. Do not immerse the scale in water or rinse it under running water.
- Store the scale in a cool, dry place when not in use; ensure no objects are placed on the scale during storage. If storing the scale for an extended period, remove the batteries.

Disposal

- Dispose of the device in accordance with local regulations. If you have any questions, please contact the local authorities responsible for waste disposal.

NOTICE

- Never use harsh solvents or cleaning products!
- Do not submerge the baby scale in water.
- Do not rinse the baby scale under running water.
- Do not clean the baby scale in the dish washer.

9. Technical Specifications

Battery type:	2x 1.5 V AA batteries
Measuring range:	0.11 lb to 44 lb (0.05 kg to 20 kg)
Weight unit:	0.005 kg , 0.01 lbs, 0.17 oz
Operating Conditions:	Temperature Range: 5°C to 35°C (41°F to 95°F); RH 20% to 95%
Storage Conditions:	Temperature Range: -20°C to +60°C (-4°F to 140°F); 20% to 95%
Scale size	22 in x 12.99 in x 2 in (559 mm x 330 mm x 51 mm)
Net weight of scale:	approx. 4.4 lb (2 kg)

Subject to technical modifications.

10. FCC Compliance Information

Baby Scale BY 80

Responsible Party – U.S. Contact Information:
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020
United States
1-800-536-0366
info@beurer.com

FCC Compliance Statement

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user‘s authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interferencein a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by

turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
—Reorient or relocate the receiving antenna.
—Increase the separation between the equipment and receiver.
—Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
—Consult the dealer or an experienced radio/TV technician for help.

11. Warranty

5-Year Limited Warranty for Original Purchaser

Your Beurer Baby Scale, Model BY 80, is warranted to be free from defects in materials and workmanship for a period of five (5) years from the date of purchase under normal conditions of intended use and service. This warranty extends only to the original retail purchaser and does not extend to retailers or subsequent owners. We will, at our option, repair or replace the Beurer Baby Scale, Model BY 80, without additional charge, for any part or parts covered by this written warranty. No refunds will be given. Repair or replacement is our only responsibility and your only remedy under this written warranty. If replacement parts for defective materials are not available, Beurer reserves the right to make product substitutions in lieu of repair or replacement.
For warranty service, contact our customer service department at 1-800-536-0366 or info@beurer.com with a description of the problem. If the issue is deemed to fall within the scope of the 5-Year Limited Warranty, you will be asked to mail the product (at your cost) in its original package along with proof of purchase, your name, address, and phone number. If the issue is not deemed eligible under warranty, we will provide a quotation for repair or replacement, including return shipping fees.

Third-Party Repairs

This warranty complies with the federal Magnuson–Moss Warranty Act. You may use third-party parts or services for repairs without voiding this warranty, unless the repair itself causes damage to the product. We encourage you to retain your proof of purchase and contact our customer service team with any questions.

Unauthorized Sellers

This warranty applies only to products purchased from Beurer or authorized retailers and distributors. Products purchased from unauthorized sellers — including third-party platforms or online marketplaces (such as eBay or unverified Amazon sellers) — are not covered by this warranty. This policy is in place to protect customers from counterfeit, altered, or used products that may not meet Beurer's quality standards.
If you're unsure whether your seller is authorized, please contact our support team.

Exclusions

This warranty does not cover damage caused by misuse or abuse; accident; the attachment of unauthorized accessories; alterations; improper installation; misapplication; lack of reasonable care; unauthorized repairs or modifications; improper electrical/power supply; worn batteries; normal wear and tear; power loss; dropped product; failure to follow instructions; lack of recommended maintenance; transit damage; theft; neglect; vandalism; or environmental conditions.
This warranty is void if the product is used in a commercial or business environment.
Beurer's maximum liability under this warranty is limited to the actual purchase price paid, as confirmed by proof of purchase. The warranty is effective only in the country where the product was purchased and used.

Legal Notice

This warranty gives you specific legal rights, and you may have other rights which vary from state to state.
THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. ANY IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ARE LIMITED IN DURATION TO THE LENGTH OF THIS WRITTEN WARRANTY.
IN NO EVENT SHALL BEURER BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING FROM USE OR INABILITY TO USE THE PRODUCT.
Beurer does not authorize any third party — including retailers or resellers — to modify the terms of this warranty. This warranty does not apply to opened, used, repackaged, or resold products, including those sold through online auctions or unauthorized channels.

For more information regarding our product line in the USA, please visit: www.shop-beurer.com

Distributed by:
Beurer North America LP
1 Oakwood Boulevard, Suite 255
Hollywood, FL 33020
United States
www.shop-beurer.com

Questions or Comments?
Call our US-based customer service toll free at 1-800-536-0366.

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